


The tumor microenvironment (TME) is comprised of a multitude of cell types that collectively create an immunosuppressive environment, enabling tumor growth. Understanding the spatial organization of these tissues requires a technology that can identify the presence of multiple markers and correlate it to their specific spatial location within the same tissue. Characteristics of the TME including infiltrating immune cells, angiogenesis and the presence of other non-malignant cells influence drug delivery, treatment effectiveness and, ultimately, clinical outcomes. Current methodologies for analyzing the spatial dimension of traditional immunofluorescence (IF) and tissues, e.g. immunohistochemistry (IHC), are limited to measuring a few parameters simultaneously, thereby restricting the number of phenotypes that can be identified. Conversely, single-cell technologies like mass cytometry and Next Generation Sequencing-based tools provide multiplexing capabilities, but at the expense of the associated spatial information.

Akoya Biosciences, Inc. is commercializing CODEX[™] (CO-Detection by indEXing), a multi-parametric imaging platform that allows the simultaneous detection and quantification of dozens of target epitopes in single cells within a single tissue section. This innovative platform is an end-to-end solution that comprises three components: 1) the CODEX fluidics instrument, 2) a suite of specialized CODEX reagents and 3) an analysis pipeline. The CODEX workflow involves a barcoding system such that each antibody moiety is conjugated to a proprietary tag (the Barcode). Panels of CODEX antibodies are used to stain tissue specimens en masse in a single step. Staining data for sets of antibodies are revealed across iterative cycles using corresponding dye-labeled Reporters. The CODEX fluidics instrument integrates with existing microscope units for a fully automated data collection process. CODEX data is processed to achieve noise reduction and analyzed at the single-cell level through a segmentation algorithm based on nuclear staining.

More than 80 antibodies have been validated on the CODEX platform for analysis of human and mouse fresh-frozen and FFPE tissues. Preliminary studies on fresh-frozen tissues with developed antibodies demonstrate an unprecedented capability for revealing the spatial correlation between a variety of target epitopes. These results demonstrate the enormous potential of the CODEX technology to identify spatial correlations in the TME through highly multiplexed detection with single-cell resolution.





A. Antibodies are conjugated to a CODEX Barcode with a highly specific corresponding dye-labeled CODEX Reporter. B. Antibody cocktail made from full panel typically of 20-50 Barcoded antibodies. C. CODEX Assay begins with single staining step with antibody cocktail (1.). Subset of targets is revealed by automatic addition of first cycle of Reporters (2.). Imaging of 3 spectrally distinct dyes for first cycle (3). Automatic removal of Cycle 1 Reporters (4.). Repeat process to assay full panel (5.) **D.** In this manner, high quality data is captured for all cycles and targets.

CODEX[®]: a novel platform for spatially-resolved deep biomarker profiling of single cells in tissue samples


Maria Elena Gallina¹, Nadya Nikulina¹, Atri Choksi¹, Jaskirat Singh¹, Gajalakshmi Dakshinamoorthy¹, Joseph Kim¹, Sejal Mistry¹, Julia Kennedy-Darling¹

¹Akoya Biosciences, Menlo Park, CA







A. The CODEX Instrument integrates with standard inverted fluorescent microscopes for automated reagent dispensation and imaging through




Clustering algorithm Xshift [1] assigned segmented cells to 16 clusters on the basis of expression levels of 21 different biomarkers. We were able to assign 15 clusters to known phenotypes present in tonsils.



The spatial location of each cell belonging to a specific phenotype (or cluster) is known, making possible to perform proximity studies and spatial analysis.









Data Clustering on CODEX data

From Clustering to Deep phenotyping in Activated T_{FH} ce Population Mature B cells Helper T cells Stromal cells •• Plasma ce Th-1 cells T_{FH} cells Plasma cells Activated T_{FH} cells DAPI nuclear staining DAPI nuclear staining

Germinal Centers





Information derived from CODEX analysis allow to build Interaction Maps, which illustrate the extent of microscopic interactions between phenotypes Here, interacting cells are defined as cells close to each other 7 μ m or less.

Conclusion

- information
- ✓ CODEX data can be analyzed through classical flow cytometry analysis platforms and/or by clustering algorithms. Data Clustering allows to identify correlation between multiple markers that are difficult to discern manually in high dimensional dataset. Clusters are manually classified as phenotypes, which can be investigated on the basis of their spatial location in the tissues, allowing us to obtain information on:
- ✓ phenotype distribution in space ✓ Microscopic interactions between single cells

References:



✓ CODEX allows to obtain high throughput single cell information coupled with Virtually Unlimited Multiplexing and Spatial

[1] N. Samusik *et al.*, Nature Methods, **2016**, 13, 493-6.